DESIGNING GOOD LIFE

# Massage REVITALIZE

*A genuine massage experience*

*Making good!*

<b>Bruksanvisning - svenska .....</b>	<b>sida 3 - 6</b>
<b>Brugsanvisning - dansk .....</b>	<b>side 7 - 10</b>
<b>Bruksanvisning - norsk .....</b>	<b>side 11 - 14</b>
<b>Käyttöohjeet - suomi.....</b>	<b>sivu 15 - 18</b>
<b>Instructions of use - english.....</b>	<b>page 19 - 22</b>

## Säkerhetsanvisningar

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## Safety instructions

1. Read the instruction manual carefully before use, and save it for future use.
2. Plug the appliance to a household outlet of 230V only, and use the appliance only for the intended use.
3. Always unplug the appliance from the mains supply when disorder is experienced, when not in use, when parts are applied or removed from the appliance and before cleaning and maintenance.
4. The appliance must not be immersed into water or other liquids.
5. This appliance can be used by children aged 8 years and above and persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, if they have been given supervision or instruction concerning the use of the appliance in a safe way by a person responsible for their safety and they understand the hazards involved.
6. Children should be supervised to ensure that they do not play with the appliance. Children are not always able to understand potential risks. Teach children a responsible handling of electrical appliances.
7. Cleaning and user maintenance shall not be made by children, unless they are over the age of 8 and supervised.
8. Check that the cord is not hanging over the edge of the table.
9. The appliance must not be covered when in use, and should only be used under supervision.

- 
- 
- 
10. Never use the massage apparatus in the bathtub, in the shower, or over a sink containing water.
  11. When in use the surface of the appliance will become warm. Persons insensitive to heat should be careful when using the appliance.
  12. Always examine appliance, cord and plug for damage before use. If the supply cord is damaged, it must be replaced by the manufacturer, or his service agent or a similarly qualified person in order to avoid a hazard.
  13. For additional protection the installation of a residual current device (RCD) with a rated residual operating current not exceeding 30 mA is recommended. Ask your installer for advice.
  14. The appliance is for domestic use only, and must not be used professionally or medically, since the appliance is not included in the CE directive concerning electro-medical appliances.
  15. If the appliance is used for other purposes than intended, or is it not handled in accordance with the instruction manual, the full responsibility for any consequences will rest with the user. Any damages to the product or other things are not covered by the warranty.

## Use

Get comfortable near a power supply. The appliance has four replaceable attachments depending on which kind of massage you want. It is recommended that you try every attachment to choose which suits best for the occasion. Screw the massage attachments on to the massage heads before you start the appliance.

Start the appliance by pressing the switch underneath the handle. The appliance always start at the highest speed (5) and in vibrating mode. Adjust the speed if you wish with the SPEED button on the handle, the lamp to each setting lights when activated. The massager gives an efficient massage and

it is recommended to start at the lowest speed (1) and increase the speed if needed.

When you have adjusted the speed you can select a program by pressing the MODE button on the handle. Selected program is indicated by respective pilot lamp. To shut off program mode, press any button. You can choose between the following program:

**BURST:** this program works in cycles with fast speed for 3 seconds and thereafter stop for 2 seconds. This will be repeated until you press any button.

**WAVE:** the speed is gradually increased from speed 2 to 5 in repeated cycles until you press any button.

**SWING:** the speed is gradually increased from speed 1 - 2 - 3 - 4 and back 4 - 3 - 2 - 1. This will be repeated until you press any button.

### Heat function

Press the HEAT button to activate the heat function. If the heat function has been activated for more than 6 minutes, make sure that the massager is not too warm. The temperature can be adjusted by shutting off the heat function for shorter periods. To shut off the heating function, press the HEAT button again.

Hold the appliance and move the appliance in circular moves over the targeted area. Avoid to press the massage heads towards tender areas.

**Back:** Move the appliance in circular moves over the back and shoulders.

**Neck:** Move the appliance in circular moves along the neck and spine.

Arms and legs: Massage in soft moves up and down along the muscles.

**Generally:** Move the appliance in circular moves to soften the muscles and for relaxation.

Do not massage the same area for more than 3 minutes, move the appliance to the next area.

Always massage in direction towards the heart no matter which area that is treated. The appliance should not be used continuously for more than 20 minutes. After 20 minutes use we recommend at least 20 minutes rest.

**Important! If you have blood circulation problems or reduced sensibility, for example if you have diabetes or neck injury, you should not use the massager before consulting your doctor. Do not use the appliance on sore, swollen or burnt skin.**

The massage should be a positive and comfortable experience and should be stopped if you feel any discomfort.

## Cleaning and maintenance

Always shut off the appliance and pull the plug before cleaning. The appliance can be wiped with a damp cloth. The appliance and its cord must not be immersed in water or other liquids.

## Disposal of the appliance

Legislation requires that all electrical and electronic equipment must be collected for reuse and recycling. Electrical and electronic equipment marked with the symbol indicating separate collection of such equipment must be returned to a municipal waste collection point.

## Warranty

Claims according to current law. Cash receipt with purchase date shall be enclosed with any claims. When raising a claim the appliance must be returned to the store of purchase.

### OBH Nordica AB

Löfströms Allé 5  
172 66 Sundbyberg  
Tel: +46 08-629 25 00  
[www.obhnordica.se](http://www.obhnordica.se)

## Technical data

OBH Nordica 6081  
230 volt  
22 watt  
Uptime: 20 min per treatment.  
4 massage attachments

With reservation for running changes.

SE/CS/6081/XX0215